

Product datasheet

PE Rat IgG2c, kappa monoclonal [SB68b] - Isotype control ab154450

1 References

Overview

Product name	PE Rat IgG2c, kappa monoclonal [SB68b] - Isotype control
Conjugation	PE. Ex: 488nm, Em: 575nm
Specificity	Keyhole limpet hemocyanin (KLH)
Tested applications	Suitable for: ELISA, Flow Cyt

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: PBS, Sucrose
Purity	Affinity purified
Clonality	Monoclonal
Clone number	SB68b
Isotype	IgG2c
Light chain type	kappa

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of **ab154450** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent concentration.
Flow Cyt		Use 0.1µg for 10 ⁶ cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors