

Product datasheet

Poly (4:1 Glu, Tyr) peptide ab204877

1 References

Overview

Product name Poly (4:1 Glu, Tyr) peptide

Description

Nature Synthetic

Specifications

Our [Abpromise guarantee](#) covers the use of **ab204877** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Functional Studies

Form Liquid

Additional notes ab204877 (Poly (4:1 Glu, Tyr) peptide) can be utilized as a substrate for the following active protein kinases:

- [ab89747](#) (Active human ACK1 protein fragment)
- [ab60876](#) (Active human BLK full length protein)
- [ab60877](#) (Active human BMX full length protein)
- [ab60764](#) (Active human CSK full length protein)
- [ab89746](#) (Active human EGFR protein fragment)
- [ab55900](#) (Active Eph receptor A1 protein fragment)
- [ab101778](#) (Active human Eph receptor A1 protein fragment)
- [ab84576](#) (Active human Eph receptor A3 protein fragment)
- [ab84575](#) (Active human Eph receptor A4 protein fragment)
- [ab56257](#) (Active mouse Eph receptor B1 protein fragment)
- [ab60865](#) (Active human Hck protein fragment)
- [ab204136](#) (Active human MATK full length protein)
- [ab72490](#) (Active human MERTK protein fragment)
- [ab60755](#) (Active PDGF Receptor alpha protein fragment)
- [ab84797](#) (Active PDGF Receptor alpha protein fragment)
- [ab60833](#) (Active PDGF Receptor beta protein fragment)
- [ab177257](#) (Active human NKF3 kinase family member protein fragment)
- [ab60878](#) (Active human PYK2 protein fragment)
- [ab60888](#) (Active human Brk full length protein)

ab107697 (Active human SRMS protein fragment)
ab105197 (Active human Tec full length protein)
ab125539 (Active human TYK2 protein fragment)
ab60334 (Active Met (c-Met) protein fragment)
ab60887 (Active TrkA protein fragment)
ab56652 (Active human TrkB protein fragment)
ab85274 (Active human TrkC protein fragment)
ab60886 (Active human Syk full length protein)
ab64304 (Active human TYRO3 protein fragment)
ab101504 (Active human Yes1 full length protein)
ab60963 (Active human ZAP70 full length protein)
ab56259 (Active human MCSF Receptor protein fragment)
ab60851 (Active Eph receptor B2 protein fragment)
ab51409 (Active human Eph receptor B3 protein fragment)
ab85713 (Active human Eph receptor B4 protein fragment)
ab126932 (Active human Eph receptor B4 protein fragment)
ab60866 (Active ErbB 2 protein fragment)
ab85602 (Active human ErbB 4 protein fragment)
ab51411 (Active human FES full length protein)
ab60853 (Active human FGFR1 protein fragment)
ab89755 (Active human FGFR2 protein fragment)
ab60857 (Active FGFR3 protein fragment)
ab89754 (Active human FGFR4 protein fragment)
ab60854 (Active human FGR full length protein)
ab126923 (Active human Flt4 protein fragment)
ab60861 (Active FRK protein fragment)
ab84696 (Active human Fyn full length protein)
ab56258 (Active mouse FER protein fragment)
ab89756 (Active human JAK3 protein fragment)
ab60305 (Active human Lck full length protein)
ab70789 (Active human Lyn full length protein)
ab84806 (Active human Met (c-Met) protein fragment)
ab105909 (Active human FAK full length protein)
ab60852 (Active human FAK protein fragment)

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors