


Product datasheet

Anti-PP4R1 antibody ab126933

1 References 1 Image

Overview

<b>Product name</b>	Anti-PP4R1 antibody
<b>Description</b>	Rabbit polyclonal to PP4R1
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>Species reactivity</b>	<b>Reacts with:</b> Human <b>Predicted to work with:</b> Mouse, Rat, Xenopus laevis, Zebrafish 
<b>Immunogen</b>	Recombinant fragment, corresponding to a region within amino acids 97-291 of of the Human PP4R1 protein (AAH60829; Q8TF05).
<b>Positive control</b>	U-87 MG whole cell lysate.
<b>General notes</b>	Keep as concentrated solution.

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.00 Preservative: 0.025% Proclin Constituents: PBS, 1% BSA, 20% Glycerol
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab126933** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/3000. Predicted molecular weight: 105 kDa.


## Target

---

<b>Relevance</b>	PP4R1 is a protein which is associated with the catalytic subunit of Protein Serine/Threonine Phosphatase 4 (PP4C) and may possibly regulate the activity of the enzyme.
<b>Cellular localization</b>	Nuclear and cytoplasmic.

## Images

---



Anti-PP4R1 antibody (ab126933) at 1/1000 dilution + U-87 MG whole cell lysate at 30 µg

Developed using the ECL technique.

**Predicted band size:** 105 kDa

7.5% SDS PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors