


Product datasheet

MG-132, proteasome inhibitor ab141003

★★★★★ [1 Abreviews](#) [32 References](#) [11 Images](#)

Overview

Product name	MG-132, proteasome inhibitor
Description	Potent, reversible proteasome inhibitor
Biological description	Potent, reversible, proteasome inhibitor ($K_i = 4 \text{ nM}$). Inhibits NF- κ B activation by preventing I κ B degradation ($IC_{50} = 3 \text{ }\mu\text{M}$). Anti-cancer properties <i>in vitro</i> and <i>in vivo</i> . Cell-permeable. Also available as ethanol solution (ab147047).
Purity	> 98%
CAS Number	133407-82-6
Chemical structure	


Properties

Molecular weight	475.62
Molecular formula	$C_{26}H_{41}N_3O_5$
Sequence	LLL (Modifications: N-terminal benzyloxycarbonyl; C-terminal aldehyde)
PubChem identifier	462382
Storage instructions	Store at -20°C . It is important to note that this is air sensitive and impurities can occur as a result of air oxidation. Store under desiccating conditions.
Solubility overview	Soluble in DMSO to 100 mM but unstable for prolonged periods. Soluble in ethanol to 100 mM.
Handling	Wherever possible, you should prepare and use solutions on the same day. However, if you need to make up stock solutions in advance, we recommend that you store the solution as aliquots in tightly sealed vials at -20°C . Generally, these will be useable for up to one week. Before use, and prior to opening the vial we recommend that you allow your product to equilibrate to room temperature for at least 1 hour. Need more advice on solubility, usage and handling? Please visit our frequently asked questions (FAQ) page for more details.
SMILES	<chem>CC(C)CC(C=O)NC(=O)C(CC(C)C)NC(=O)C(CC(C)C)NC(=O)OCC1=CC=CC=C1</chem>
Source	Synthetic



Chemical Structure - MG-132, proteasome inhibitor
(ab141003)

2D chemical structure image of ab141003, MG-132, proteasome inhibitor



Western blot - MG-132, proteasome inhibitor
(ab141003)

All lanes : Anti-SQSTM1 / p62 (phospho S349) antibody [EPR20451] (**ab211324**) at 1/1000 dilution

Lane 1 : Untreated HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate

Lane 2 : HeLa whole cell lysate treated with 2 μM MG-132 (ab141003) for 18 hours

Lane 3 : HeLa whole cell lysate treated with 2 μM MG-132 (ab141003) for 18 hours, then treated with Alkaline Phosphatase for 1 hour

Lysates/proteins at 10 μg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/100000 dilution

Observed band size: 62 kDa

Exposure time: 1 minute

Blocking/Dilution buffer: 5% NFDN/TBST.



Immunocytochemistry/ Immunofluorescence - MG-132, proteasome inhibitor (ab141003)


Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized HeLa (Human epithelial cell line from cervix adenocarcinoma) cells, treated with 2µM MG-132 (ab141003) for 18 hours or untreated, labeling SQSTM1 / p62 (phospho S349) with **ab211324** at 1/100 dilution, followed by Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) (**ab195889**) secondary antibody at 1/200 dilution (green).

Confocal image showing cytoplasmic staining on HeLa cell line. The expression increased after treatment with 2µM MG-132 (ab141003) for 18 hours.

The nuclear counterstain is DAPI (blue).


Tubulin is detected with **ab195889** (Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594)) at 1/200 dilution (red).

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat anti-rabbit IgG (Alexa Fluor® 488) (**ab150077**) at 1/1000 dilution.




Immunocytochemistry/ Immunofluorescence - MG-132, proteasome inhibitor (ab141003)

ab62352 staining Nrf2 in untreated HeLa cells (top panel) and treated HeLa cells (bottom panel). Cells were treated with 2µM of MG-132 for 18 hours (ab141003). The cells were fixed with 4% paraformaldehyde (10 min) then permeabilized with 0.1% Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated with **ab62352** at 1µg/ml concentration and **ab7291** (Mouse monoclonal to alpha Tubulin) at 1/1000 dilution overnight at 4°C followed by a further incubation at room temperature for 1h with a goat secondary antibody to rabbit IgG (Alexa Fluor® 488) (**ab150081**) at 2µg/ml (shown in green) and a goat secondary antibody to mouse IgG (Alexa Fluor® 594) (**ab150120**) at 2µg/ml (shown in red). Nuclear DNA was labelled in blue with DAPI. Image was taken with a confocal microscope (Leica-Microsystems TCS SP8).



Flow Cytometry (Intracellular) - MG-132, proteasome inhibitor (ab141003)

Flow cytometric analysis of 4% paraformaldehyde-fixed HeLa (Human epithelial cell line from cervix adenocarcinoma) cells, treated with 2µM MG-132 (ab141003) for 18 hours (red) or untreated (green), labeling SQSTM1 / p62 (phospho S349) with **ab211324** at 1/500 dilution compared with a rabbit monoclonal IgG isotype control (**ab172730**; black) and an unlabeled control (cells without incubation with primary antibody and secondary antibody; blue). Goat anti rabbit IgG (Alexa Fluor® 488) at 1/2000 dilution was used as the secondary antibody.



Immunoprecipitation - MG-132, proteasome inhibitor (ab141003)

SQSTM1 / p62 (phospho S349) was immunoprecipitated from 0.35 mg of HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate treated with 2µM MG-132 (ab141003) for 18h with **ab211324** at 1/30 dilution.

Western blot was performed from the immunoprecipitate using **ab211324** at 1/1000 dilution.


VeriBlot for IP Detection Reagent (HRP) (**ab131366**), was used for detection at 1/10000 dilution.

Lane 1: HeLa treated with 2µM MG-132 (ab141003) for 18h whole cell lysate, 10 µg (Input).

Lane 2: **ab211324** IP in HeLa treated with 2µM MG-132 (ab141003) for 18h whole cell lysate.

Lane 3: Rabbit monoclonal IgG (**ab172730**) instead of **ab211324** in HeLa treated with 2µM MG-132 (ab141003) for 18h whole cell lysate.

Blocking and dilution buffer: 5% NFDM/TBST.



All lanes : Anti-SQSTM1 / p62 (phospho S349) antibody [EPR20451] ([ab211324](#)) at 1/1000 dilution

Lane 1 : Untreated C6 (Rat glial tumor cell line) whole cell lysate

Lane 2 : C6 whole cell lysate treated with 2 μ M MG-132 (ab141003) for 18 hours

Lysates/proteins at 10 μ g per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/100000 dilution

Observed band size: 62 kDa

Exposure time: 1 minute

Blocking/Dilution buffer: 5% NFDN/TBST.



All lanes : Anti-SQSTM1 / p62 (phospho S349) antibody [EPR20451] ([ab211324](#)) at 1/1000 dilution

Lane 1 : Untreated NIH/3T3 (Mouse embryonic fibroblast cell line) whole cell lysate

Lane 2 : NIH/3T3 whole cell lysate treated with 2 μ M MG-132 (ab141003) for 18 hours

Lysates/proteins at 10 μ g per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/100000 dilution

Observed band size: 62 kDa

Exposure time: 4 seconds

Blocking/Dilution buffer: 5% NFDM/TBST.

All lanes : Anti-HIF-1 alpha antibody [EPR16897] (**ab179483**) at 0.163 µg/ml

Lane 1 : Untreated C6 (rat glial tumor glial cell), whole cell lysate

Lane 2 : C6 treated with 400 µM CoCl₂ and 20 µM MG-132 (ab141003) for 4 hours

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

Observed band size: 110 kDa

Exposure time: 26 seconds

Blocking and diluting buffer: 5% NFDM/TBST.

The expression of HIF-1 alpha is induced by CoCl₂ and maintained by MG-132 (PMID: 15836611).

All lanes : Anti-YTHDF2 antibody [EPR20318] (**ab220163**) at 1/5000 dilution

Lane 1 : Untreated HT-1080 (human fibrosarcoma epithelial cell) whole cell lysate at 20 µg


Lane 2 : HT-1080 treated with 10µM MG-132 (ab141003) for 4 hours, whole cell lysate at 10 µg

Secondary


All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/100000 dilution

Exposure time: 103 seconds


Blocking/diluting buffer and concentration: 5% NFDM/TBST.



Western blot - MG-132, proteasome inhibitor (ab141003)



Western blot - MG-132, proteasome inhibitor (ab141003)



Western blot - MG-132, proteasome inhibitor (ab141003)

All lanes : Anti-ATF-4 antibody [EPR18111] (**ab184909**) at 1/1000 dilution

Lane 1 : Untreated HepG2 (Human liver hepatocellular carcinoma cell line) whole cell lysate (control)

Lane 2 : HepG2 (Human liver hepatocellular carcinoma cell line) treated with 5 μ M MG-132 (ab141003) and 3 μ g/ml tunicamycin (**ab120296**) for 6 hours whole cell lysate

Lysates/proteins at 20 μ g per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/100000 dilution

Observed band size: 50 kDa

Exposure time: 15 seconds

Blocking/Dilution buffer: 5% NFDm/TBST.

The molecular weight observed is consistent with what has been described in the literature (PMID: 22095285).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES, NOT FOR USE IN HUMANS"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors
- Abcam biochemicals are novel compounds and we have not tested their biological activity in house. Please use the literature to

identify how to use these products effectively. If you require further assistance please contact the scientific support team