

Product datasheet

Pluronic® F-127 ab285413

Overview

Product name	Pluronic® F-127
Description	Non-ionic detergent for isolation of membrane proteins
General notes	This product is manufactured by BioVision, an Abcam company and was previously called 2730 Pluronic® F-127. 2730-250G is the same size as the 250 g size of ab285413.
CAS Number	9003-11-6

Properties

Chemical name	2-methyloxirane;oxirane
Molecular weight	12500.00
Molecular formula	C ₅ H ₁₀ O ₂
PubChem identifier	24751
Storage instructions	Shipped at Room Temperature. Store at Room Temperature. Store under desiccating conditions.
Solubility overview	H ₂ O (50 mg/ml)
Handling	Need more advice on solubility, usage and handling? Please visit our frequently asked questions (FAQ) page for more details.
SMILES	CC1CO1.C1CO1

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES, NOT FOR USE IN HUMANS"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors
- Abcam biochemicals are novel compounds and we have not tested their biological activity in house. Please use the literature to identify how to use these products effectively. If you require further assistance please contact the scientific support team