


Product datasheet

Polyethylene Glycol Dimethacrylate ab254480

Overview

Product name	Polyethylene Glycol Dimethacrylate
Description	Polymer used in hydrogels.
Biological description	Polymer used in a range of chemical processes. Used in hydrogels. Used as adhesives, coatings, sealants, photo resists, solder masks and photopolymers.
CAS Number	25852-47-5
Chemical structure	

Properties

Storage instructions	Shipped at Room Temperature. Store at Room Temperature.
Handling	Refer to SDS for further information. Need more advice on solubility, usage and handling? Please visit our frequently asked questions (FAQ) page for more details.
Source	Synthetic

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES, NOT FOR USE IN HUMANS"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors
- Abcam biochemicals are novel compounds and we have not tested their biological activity in house. Please use the literature to identify how to use these products effectively. If you require further assistance please contact the scientific support team