

Product datasheet

Normal Rat Serum (sterile) ab7488

4 References

Overview

Product name	Normal Rat Serum (sterile)
Host species	Rat
Tested applications	Suitable for: Blocking, IHC-P, IHC-Fr
General notes	<p>Rat serum for blocking in antibody-based applications or use as an antibody diluent. From mixed sex Sprague Dawley rats. Origin: North America.</p> <p>Total Protein Concentration : normal.</p> <p>Purity: Whole antiserum</p>

Properties

Form	Liquid
Storage instructions	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab7488 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Blocking		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors