

Human wild-type SH-SY5Y cell line **ab275475**

1 References

Overview

Product name	Human wild-type SH-SY5Y cell line
Parental Cell Line	SHSY-5Y
Organism	Human
Passage number	<20
Biosafety level	1
General notes	Wild-type cell lines are sold with knockout cell lines only - not available for individual purchase.

Cryopreservation cell medium: Cell Freezing Medium-DMSO Serum free media, contains 8.7% DMSO in MEM supplemented with methyl cellulose.

Culture medium: 1:1 mixture of EMEM and F-12K + 10% FBS

Initial handling guidelines: Upon arrival, the vial should be stored in liquid nitrogen vapor phase and not at -80°C. Storage at -80°C may result in loss of viability.

1. Thaw the vial in 37°C water bath for approximately 1-2 minutes.
2. Transfer the cell suspension (0.8 mL) to a 15 mL/50 mL conical sterile polypropylene centrifuge tube containing 8.4 mL pre-warmed culture medium, wash vial with an additional 0.8 mL culture medium (total volume 10 mL) to collect remaining cells, and centrifuge at 201 x g (rcf) for 5 minutes at room temperature. 10 mL represents minimum recommended dilution. 20 mL represents maximum recommended dilution.
3. Resuspend the cell pellet in 5 mL pre-warmed culture medium and count using a haemocytometer or alternative cell counting method. Based on cell count, seed cells in an appropriate cell culture flask at a density of $1-2 \times 10^4$ cells/cm². Seeding density is given as a guide only and should be scaled to align with individual lab schedules.
4. Incubate the culture at 37°C incubator with 5% CO₂. Cultures should be monitored daily.

Subculture guidelines:

All seeding densities should be based on cell counts gained by established methods.

These cells grow as a mixture of floating and adherent cells.

Remove media containing floating cells and recover cells by centrifugation, detach cells using standard methods, combine with floating cells and transfer to a new culture flask.

A guide seeding density of $1-2 \times 10^4$ cells/cm² is recommended.

Cells should be seeded at a density conducive to cell-cell communication to proliferate. If cells are seeded too sparsely, growth rate is reduced and cell death is high.

We will provide viable cells that proliferate on revival.

Properties

Number of cells	1 x 10 ⁶ cells/vial, 1 mL
Adherent /Suspension	Adherent
Tissue	Bone marrow
Cell type	neuroblastoma
Disease	Neuroblastoma
Gender	Female
Mycoplasma free	Yes
Storage instructions	Shipped on Dry Ice. Store in liquid nitrogen.
Storage buffer	Constituents: 8.7% Dimethylsulfoxide, 2% Cellulose, methyl ether

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors