

Product datasheet

StayGreen/AP Plus (Alcohol and Xylene Substitute Compatible) ab156428

★★★★★ 1 Abreviews 1 References 2 Images

Overview


Product name	StayGreen/AP Plus (Alcohol and Xylene Substitute Compatible)
Assay time	0h 20m
Product overview	<p>StayGreen/AP is a substrate chromogen system designed to be used for either IHC or ISH when using alkaline phosphate detection. StayGreen/AP produces a green color. It is insoluble in alcohol and xylene solutions (both aliphatic hydrocarbon and citrus based) therefore sections can be dehydrated in alcohol, cleared in xylene substitute and permanently mounted.</p> <p>This chromogen substrate system may be used for both automation and manual use.</p>

Properties

Storage instructionsStore at +4°C. Please refer to protocols.


Components	1 kit
Clear StayGreen/AP Substrate Buffer	1 x 30ml
Concentrated StayGreen/ AP Chromogen	1 x 1ml
Empty Mixing Dropping Bottle	1 unit

Images



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - StayGreen/AP (Alcohol and Xylene Substitute Compatible) (ab156428)

Formalin fixed paraffin embedded Human tonsil stained with HMW Cytokeratin antibody ([ab850](#)) labeled with StayGreen/AP (Alcohol and Xylene Substitute Compatible) (ab156428) produces a strong green color.



Example of Staining with ab156428

Formalin-fixed, paraffin-embedded human tonsil stained with HMW CK antibody labeled with StayGreen/AP (ab156428) produces a strong green color.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors