

Anti-ANGPTL3 antibody [EPR21979-129] - BSA and Azide free (Capture/Detector) ab259731

Recombinant RabMAb

3 Images

Overview

Product name	Anti-ANGPTL3 antibody [EPR21979-129] - BSA and Azide free (Capture/Detector)
Description	Rabbit monoclonal [EPR21979-129] to ANGPTL3 - BSA and Azide free (Capture/Detector)
Host species	Rabbit
Tested applications	Suitable for: Sandwich ELISA
Species reactivity	Reacts with: Mouse, Human
Immunogen	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers.
General notes	<p>ab259731 is a BSA and Azide Free antibody supplied in an unconjugated format and it is suitable for sandwich ELISAs to quantify Mouse ANGPTL3 (Angiopietin-like 3). The recommended pair for sandwich ELISA is:</p> <p>Capture: ab281289, Mouse ANGPTL3 (Angiopietin-like 3) Capture Antibody (unconjugated)</p> <p>Detector: ab259731, Mouse ANGPTL3 (Angiopietin-like 3) Detector Antibody (unconjugated)</p> <p>The reference range value is 125-4000 pg/mL.</p> <p>Alternatively, you can use:</p> <p>Capture: ab259731, Human ANGPTL3 (Angiopietin-like 3) Capture Antibody (unconjugated)</p> <p>Detector: ab259732, Human ANGPTL3 (Angiopietin-like 3) Detector Antibody (unconjugated)</p> <p>The reference range value is 218.75 - 14000 pg/ml.</p> <p>Our carrier-free antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency.</p> <p>This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications.</p> <p>Use our conjugation kits for antibody conjugates that are ready-to-use in as little as 20 minutes with <1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.</p> <p>The recommended antibody orientation is based on internal optimization for ELISA-based assays. Antibody orientation is assay dependent and needs to be optimized for each assay type. Please note that the range provided for this antibody is only an estimation based on the</p>

performance of the product using the recommended antibody pair. Performance of the antibody pair will depend on the specific characteristics of your assay. We guarantee the product works in sandwich ELISA, but we do not guarantee the sensitivity or dynamic range of the antibody in your assay.

This product is a recombinant monoclonal antibody, which offers several advantages including:

- High batch-to-batch consistency and reproducibility
- Improved sensitivity and specificity
- Long-term security of supply
- Animal-free production

For more information [see here](#).

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb[®] patents](#).

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Constituent: 100% PBS
Carrier free	Yes
Purity	Protein A purified
Clonality	Monoclonal
Clone number	EPR21979-129
Isotype	IgG

Applications


The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab259731 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Sandwich ELISA		Use at an assay dependent concentration. Can be paired for Sandwich ELISA with Rabbit monoclonal [EPR21979-110] to ANGPTL3 - BSA and Azide free (Detector) (ab259732) and Rabbit monoclonal [EPR21979-47] to ANGPTL3 - BSA and Azide free (Capture) (ab281289) .

Target


Tissue specificity	Expressed principally in liver. Weakly expressed in kidney.
Sequence similarities	Contains 1 fibrinogen C-terminal domain.
Cellular localization	Secreted.

Images



Sandwich ELISA - Anti-ANGPTL3 antibody
[EPR21979-129] - BSA and Azide free
(Capture/Detector) (ab259731)

Sandwich ELISA of **ab253729** with the capture antibody **ab281289** dilution at 2 µg/mL and detector antibody ab259731 dilution at 0.5 µg/mL.



Sandwich ELISA - Anti-ANGPTL3 antibody
[EPR21979-129] - BSA and Azide free
(Capture/Detector) (ab259731)

Representative standard curve from corresponding SimpleStep ELISA® Kit (**ab254510**).

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-ANGPTL3 antibody [EPR21979-129] - BSA and Azide free (Capture/Detector) (ab259731)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors