


## Product datasheet

### Anti-APC2 antibody ab113370

★★★★★ [1 Abreviews](#) [3 References](#) [3 Images](#)

#### Overview

<b>Product name</b>	Anti-APC2 antibody
<b>Description</b>	Rabbit polyclonal to APC2
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> WB, IHC-P
<b>Species reactivity</b>	<b>Reacts with:</b> Human <b>Predicted to work with:</b> Mouse 
<b>Immunogen</b>	Synthetic peptide corresponding to Human APC2 (internal sequence) conjugated to keyhole limpet haemocyanin. Database link: <a href="#">O95996</a>
<b>Positive control</b>	Human Brain Cortex and Small Intestine tissues; HeLa cell lysate.
<b>General notes</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.4 Preservative: 0.02% Sodium azide Constituent: 99% PBS
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

## Applications

### The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab113370 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 244 kDa.
IHC-P	★★★★★ (1)	Use a concentration of 5 µg/ml.

## Target

### Function

Promotes rapid degradation of CTNNB1 and may function as a tumor suppressor. May function in Wnt signaling.

### Tissue specificity

Widely expressed (at protein level). Specifically expressed in the CNS.

### Sequence similarities

Belongs to the adenomatous polyposis coli (APC) family.  
Contains 6 ARM repeats.


### Developmental stage

Expressed in fetal brain.

### Cellular localization


Cytoplasm, cytoskeleton. Golgi apparatus. Cytoplasm. Cytoplasm, perinuclear region. Associated with actin filaments (PubMed:11691822). Associated with microtubule network (PubMed:10644998, PubMed:11691822).

## Images




ab113370, at 5 µg/ml, staining APC2 in Formalin-fixed, Paraffin-embedded Human Brain cortex tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-APC2 antibody (ab113370)



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-APC2 antibody (ab113370)

ab113370, at 5 µg/ml, staining APC2 in Formalin-fixed, Paraffin-embedded Human Small intestine tissue by Immunohistochemistry.



Western blot - Anti-APC2 antibody (ab113370)

**Lane 1** : Anti-APC2 antibody (ab113370) at 1 µg/ml

**Lane 2** : Anti-APC2 antibody (ab113370) at 2 µg/ml

**All lanes** : HeLa cell lysate

**Predicted band size:** 244 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors