


Product datasheet

Anti-APG8A/ATG8A antibody ab77003

★★★★★ 1 Abreviews 30 References

Overview

Product name	Anti-APG8A/ATG8A antibody
Description	Rabbit polyclonal to APG8A/ATG8A
Host species	Rabbit
Species reactivity	Reacts with: Arabidopsis thaliana, Corn Predicted to work with: Human 
Immunogen	Recombinant full length protein corresponding to APG8A/ATG8A.
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Constituents: 2.42% Tris, 1.5% Glycine
Purity	IgG fraction
Purification notes	This antibody was purified over AFFI-GEL 10 beads binds with rATG8 protein.
Clonality	Polyclonal
Isotype	IgG

Target

Relevance	APG8A is involved in cytoplasm to vacuole transport (Cvt) vesicles and autophagosomes formation. It may mediate the delivery of the vesicles and autophagosomes to the vacuole via the microtubule cytoskeleton.
------------------	--

Cellular localization

Cytoplasmic vesicle › cvt vesicle membrane; Lipid-anchor. Cytoplasmic vesicle › autophagosome membrane; Lipid-anchor. Vacuole membrane; Lipid-anchor. Note= Membrane-associated through a lipid anchor.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors