

# FITC Rat IgG1, kappa [RTK2071] - isotype control ab18404

### Overview

<b>Product name</b>	FITC Rat IgG1, kappa [RTK2071] - isotype control
<b>Conjugation</b>	FITC. Ex: 493nm, Em: 528nm
<b>Conjugation notes</b>	The antibody was conjugated with FITC under optimal conditions. The solution is free of unconjugated FITC.
<b>Specificity</b>	The RTK2071 immunoglobulin reacts with KLH. The RTK2071 immunoglobulin was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse and human tissues.
<b>Tested applications</b>	<b>Suitable for:</b> Flow Cyt
<b>Immunogen</b>	Chemical/ Small Molecule corresponding to Rat conjugated to keyhole limpet haemocyanin.
<b>General notes</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituent: PBS</p>
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	RTK2071
<b>Isotype</b>	IgG1
<b>Light chain type</b>	kappa

## Applications

---

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab18404 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors