

## Product datasheet

# Anti-FOXA1 antibody [EPR10881] - ChIP Grade ab170933

**KO VALIDATED** Recombinant RabMAB

★★★★★ [1 Abreviews](#) [28 References](#) [12 Images](#)

### Overview

<b>Product name</b>	Anti-FOXA1 antibody [EPR10881] - ChIP Grade
<b>Description</b>	Rabbit monoclonal [EPR10881] to FOXA1 - ChIP Grade
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> Flow Cyt (Intra), WB, IHC-P, ICC/IF <b>Unsuitable for:</b> ChIP
<b>Species reactivity</b>	<b>Reacts with:</b> Mouse, Rat, Human
<b>Immunogen</b>	Recombinant fragment within Human FOXA1 aa 350 to the C-terminus. The exact sequence is proprietary. Database link: <a href="#">P55317</a>
<b>Positive control</b>	WB: HeLa, Hap1, SW480 and HepG2 whole cell lysate ( <a href="#">ab7900</a> ). Mouse and rat lung lysates IHC-P: Human breast carcinoma & prostate tissue, mouse liver, Rat pancreas tissue. ICC/IF: PC-3 and HepG2 cells.
<b>General notes</b>	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"><li>- High batch-to-batch consistency and reproducibility</li><li>- Improved sensitivity and specificity</li><li>- Long-term security of supply</li><li>- Animal-free production</li></ul> <p>For more information <a href="#">see here</a>.</p> <p>Our RabMAB<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAB<sup>®</sup> patents</a>.</p> <p><b>We are constantly working hard to ensure we provide our customers with best in class antibodies. As a result of this work we are pleased to now offer this antibody in purified format. We are in the process of updating our datasheets. The purified format is designated 'PUR' on our product labels. If you have any questions regarding this update, please contact our Scientific Support team.</b></p>

### Properties

**Form** Liquid

<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Preservative: 0.01% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), 59% PBS, 0.05% BSA
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	EPR10881
<b>Isotype</b>	IgG

## Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab170933 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt (Intra)		Use at an assay dependent concentration.
WB	★★★★★ (1)	1/1000 - 1/10000. Predicted molecular weight: 49 kDa.
IHC-P		1/1000. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. See <a href="#">IHC antigen retrieval protocols</a> . <b>For unpurified use 1/100 - 1/250</b>
ICC/IF		1/100 - 1/250.

**Application notes** Is unsuitable for ChIP.

## Target

**Function** Transcription factor that is involved in embryonic development, establishment of tissue-specific gene expression and regulation of gene expression in differentiated tissues. Is thought to act as a 'pioneer' factor opening the compacted chromatin for other proteins through interactions with nucleosomal core histones and thereby replacing linker histones at target enhancer and/or promoter sites. Binds DNA with the consensus sequence 5'-[AC]A[AT]T[AG]TT[GT][AG][CT]T[CT]-3' (By similarity). Proposed to play a role in translating the epigenetic signatures into cell type-specific enhancer-driven transcriptional programs. Its differential recruitment to chromatin is dependent on distribution of histone H3 methylated at 'Lys-5' (H3K4me2) in estrogen-regulated genes. Involved in the development of multiple endoderm-derived organ systems such as liver, pancreas, lung and prostate; FOXA1 and FOXA2 seem to have at least in part redundant roles (By similarity). Modulates the transcriptional activity of nuclear hormone receptors. Is involved in ESR1-mediated transcription; required for ESR1 binding to the NKX2-1 promoter in breast cancer cells; binds to the RPRM promoter and is required for the estrogen-induced repression of RPRM. Involved in regulation of apoptosis by inhibiting the expression of BCL2. Involved in cell cycle regulation by activating expression of CDKN1B, alone or in conjunction with BRCA1. Originally described as a transcription activator for a number of liver genes such as AFP, albumin, tyrosine aminotransferase, PEPCK, etc. Interacts with the cis-acting regulatory regions of these

genes. Involved in glucose homeostasis.

### Tissue specificity

Highly expressed in prostate and ESR1-positive breast tumors. Overexpressed in esophageal and lung adenocarcinomas.


### Sequence similarities

Contains 1 fork-head DNA-binding domain.

### Cellular localization

Nucleus.

## Images



Western blot - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)

**All lanes :** Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933) at 1/1000 dilution

**Lane 1 :** Wild-type HAP1 whole cell lysate

**Lane 2 :** FOXA1 knockout HAP1 whole cell lysate

**Lane 3 :** MCF7 whole cell lysate


**Lane 4 :** HepG2 whole cell lysate

Lysates/proteins at 20 µg per lane.

**Predicted band size:** 49 kDa

**Lanes 1 - 4:** Merged signal (red and green). Green - ab170933 observed at 52 kDa. Red - loading control, **ab9484**, observed at 37 kDa.

Unpurified ab170933 was shown to specifically react with FOXA1 in wild-type HAP1 cells as signal was lost in FOXA1 knockout cells. Wild-type and FOXA1 knockout samples were subjected to SDS-PAGE. Ab170933 and **ab9484** (Mouse anti-GAPDH loading control) were incubated overnight at 4°C at 1/1000 dilution and 1/20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preabsorbed **ab216773** and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed **ab216776** secondary antibodies at 1/10000 dilution for 1 hour at room temperature before imaging.



Western blot - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)

**All lanes** : Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933) at 1/1000 dilution

**Lane 1** : Wild-type HeLa cell lysate

**Lane 2** : FOXA1 knockout HeLa cell lysate

Lysates/proteins at 20 µg per lane.


Performed under reducing conditions.

**Predicted band size:** 49 kDa

**Observed band size:** 52 kDa

**Lanes 1-2:** Merged signal (red and green). Green - ab170933 observed at 52 kDa. Red - Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) observed at 37 kDa.

ab170933 was shown to react with FOXA1 in wild-type HeLa cells in western blot. Loss of signal was observed when knockout cell line [ab261823](#) (knockout cell lysate [ab256920](#)) was used. Wild-type HeLa and FOXA1 knockout HeLa cell lysates were subjected to SDS-PAGE. Membrane was blocked for 1 hour at room temperature in 0.1% TBST with 3% non-fat dried milk. ab170933 and Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) overnight at 4°C at a 1 in 1000 dilution and a 1 in 20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye®800CW) preadsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye®680RD) preadsorbed ([ab216776](#)) secondary antibodies at 1 in 20000 dilution for 1 hour at room temperature before imaging.



**All lanes :** Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933) at 1.1 µg/ml (purified)

**Lane 1 :** SW480 (Human colorectal adenocarcinoma epithelial cell) whole cell lysates

**Lane 2 :** Mouse lung lysates

**Lane 3 :** Rat lung lysates

Lysates/proteins at 15 µg per lane.


### Secondary

**All lanes :** Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

**Predicted band size:** 49 kDa


**Observed band size:** 49 kDa

Blocking and diluting buffer: 5% NFDm/TBST.




Immunocytochemistry/ Immunofluorescence analysis of PC-3 (Human prostate adenocarcinoma epithelial cell) cells labeling FOXA1 with Purified ab170933 at 1:100 dilution (11 µg/ml). Cells were fixed in 4% Paraformaldehyde and permeabilized with 0.1% tritonX-100. Cells were counterstained with Ab195889 Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor®594) 1:200 (2.5 µg/ml). Goat anti rabbit IgG (Alexa Fluor®488, **ab150077**) was used as the secondary antibody at 1:1000 (2 µg/ml) dilution. DAPI nuclear counterstain. PBS instead of the primary antibody was used as the secondary antibody only control.

Immunocytochemistry/ Immunofluorescence - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Rat pancreas tissue sections labeling FOXA1 with Purified ab170933 at 1:1000 dilution (1.13  $\mu\text{g/ml}$ ). Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0) ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Mouse liver tissue sections labeling FOXA1 with Purified ab170933 at 1:1000 dilution (1.13  $\mu\text{g/ml}$ ). Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0) ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Human breast cancer tissue sections labeling FOXA1 with Purified ab170933 at 1:1000 dilution (1.13 µg/ml). Heat mediated antigen retrieval was performed using [ab93684](#) (Tris/EDTA buffer, pH 9.0) ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain.




Immunocytochemistry/ Immunofluorescence - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)

Immunofluorescence analysis of HepG2 cells labeling FOXA1 with unpurified ab170933 at 1/100 dilution.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)

Immunohistochemical analysis of paraffin-embedded Human breast carcinoma tissue labeling FOXA1 with unpurified ab170933 at 1/100 dilution. Heat mediated antigen retrieval was performed using citrate buffer pH 6.0 before commencing with IHC staining protocol.



Flow Cytometry (Intracellular) - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)

Intracellular Flow Cytometry analysis of HepG2 (human hepatocellular carcinoma) cells labeling FOXA1 with purified ab170933 at 1/200 dilution (10ug/ml) (red). Cells were fixed with 4% paraformaldehyde and permeabilised with 90% methanol. A Goat anti rabbit IgG (Alexa Fluor® 488) (1/2000 dilution) was used as the secondary antibody. Rabbit monoclonal IgG (Black) was used as the isotype control, cells without incubation with primary antibody and secondary antibody (Blue) were used as the unlabeled control.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)

Immunohistochemical analysis of paraffin-embedded Human prostate tissue labeling FOXA1 with unpurified ab170933 at 1/100 dilution. Heat mediated antigen retrieval was performed using citrate buffer pH 6.0 before commencing with IHC staining protocol.

**Why choose a recombinant antibody?**

 <p><b>Research with confidence</b> Consistent and reproducible results</p>	 <p><b>Long-term and scalable supply</b> Recombinant technology</p>
 <p><b>Success from the first experiment</b> Confirmed specificity</p>	 <p><b>Ethical standards compliant</b> Animal-free production</p>

Anti-FOXA1 antibody [EPR10881] - ChIP Grade (ab170933)



**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors