




Anti-GnRHR antibody [F1G4] ab24098

Overview

Product name	Anti-GnRHR antibody [F1G4]
Description	Mouse monoclonal [F1G4] to GnRHR
Host species	Mouse
Specificity	It is unknown if ab24098 reacts with the receptor form 1 or 2.
Tested applications	Suitable for: ELISA
Species reactivity	Reacts with: Human Predicted to work with: Rabbit, Pig 
Immunogen	Synthetic peptide: MANSASPEQNQHCSAINNSIPLMQGNLPY conjugated to BSA, corresponding to N terminal extracellular domain amino acids 1-29 of Human GnRHR <div>  Run BLAST with  Run BLAST with </div>
Positive control	IHC-Fr: human anterior pituitary WB/FACS: GT7-1 cells (hypothalamus cell line), T47D breast carcinoma cells
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
Storage buffer	pH: 7.3 Preservative: 0.05% Sodium azide Constituent: PBS
Purity	Protein G purified

Clonality	Monoclonal
Clone number	F1G4
Myeloma	Sp2/0
Isotype	IgG1

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab24098 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent concentration.

Target

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors