

Product datasheet

Anti-hapten 4-hydroxy-3-nitrophenyl acetyl (NP) antibody [B1-8] ab206532

Recombinant

[1 References](#) [1 Image](#)

Overview

Product name	Anti-hapten 4-hydroxy-3-nitrophenyl acetyl (NP) antibody [B1-8]
Description	Rabbit monoclonal [B1-8] to hapten 4-hydroxy-3-nitrophenyl acetyl (NP)
Host species	Rabbit
Tested applications	Suitable for: WB, Flow Cyt
Species reactivity	Reacts with: Species independent
Immunogen	Full length protein corresponding to hapten 4-hydroxy-3-nitrophenyl acetyl (NP). Hapten-protein conjugated NP12-chicken gammaglobulin (CG).
General notes	ab206532 is a recombinant monoclonal antibody to 4-hydroxy-3-nitrophenyl acetyl (NP), manufactured using variable regions (i.e. specificity) from the hybridoma B1-8. It is a chimeric antibody created as part of a panel offering antibodies of the same specificity in different formats (species, isotype, subtype and modified versions) for use as isotype controls.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.02% Proclin 300 Constituent: 99% PBS
Purity	Protein A purified
Clonality	Monoclonal
Clone number	B1-8
Isotype	IgG
Light chain type	lambda

Applications

The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab206532 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration.

Images

Why choose a recombinant antibody?

- Research with confidence**
Consistent and reproducible results
- Long-term and scalable supply**
Recombinant technology
- Success from the first experiment**
Confirmed specificity
- Ethical standards compliant**
Animal-free production

Anti-hapten 4-hydroxy-3-nitrophenyl acetyl (NP) antibody [B1-8] (ab206532)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors