

HRP Anti-polyHistidine antibody [HIS-1] ab49781

2 References

Overview

| | |
|----------------------------|--|
| Product name | HRP Anti-polyHistidine antibody [HIS-1] |
| Description | HRP Mouse monoclonal [HIS-1] to polyHistidine |
| Host species | Mouse |
| Conjugation | HRP |
| Tested applications | Suitable for: ELISA, Dot blot, WB |
| Species reactivity | Reacts with: Species independent |
| Immunogen | Recombinant polyHistidine tagged fusion protein. |
| Positive control | Bacteria lysates expressing a recombinant Histidine tagged fusion protein. |
| General notes | <p>Monoclonal antibody reacting specifically with polyHistidine may be useful in various immunotechniques, to identify the expression of a polyhistidine fusion protein in bacteria, bacterial lysates or cells and tissues transfected with a polyHistidine tagged fusion protein expressing vectors.</p> <p>The conjugate was prepared by conjugation of HRP to Mouse monoclonal anti-polyHistidine purified from ascites fluid of the HIS-1 hybridoma. Molar ratio Ab/Enzyme: 0.8 to 1.5 Enzyme activity: at least 400 U/ml</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C. |
| Storage buffer | Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.0268% PBS, 1% BSA |
| Purity | IgG fraction |

| | |
|-------------------------------|---|
| Purification notes | Purified Immunoglobulin |
| Primary antibody notes | Monoclonal antibody reacting specifically with polyHistidine may be useful in various immunotechniques, to identify the expression of a polyhistidine fusion protein in bacteria, bacterial lysates or cells and tissues transfected with a polyHistidine tagged fusion protein expressing vectors. |
| Clonality | Monoclonal |
| Clone number | HIS-1 |
| Isotype | IgG2a |

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab49781 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-----------------|-----------|--|
| ELISA | | Use at an assay dependent concentration. |
| Dot blot | | 1/2000. |
| WB | | 1/2000. |

Target

Relevance Recombinant DNA technology enables the attachment of genes of interest to specific sequences or genes that can provide 'affinity handles' (tags) designed to enable the selective identification of the protein of interest. These sequences of tails or tags are genetically engineered away from the protein active site, by insertion at the N or C terminus. It has been reported that the addition of a consecutive Histidine amino acid residue tail creates a stable fusion product that does not appear to interfere with the bioactivity of the protein or with the biodistribution of the Histidine tagged product.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors