

## Product datasheet

# Human IgG - Isotype Control ab206195

## 2 References

### Overview

**Product name** Human IgG - Isotype Control

**Tested applications** **Suitable for:** ELISA

**General notes**

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

### Properties

**Form** Liquid

**Storage instructions** Shipped at 4°C. Store at +4°C.

**Storage buffer** pH: 8.2  
Constituents: 0.6% Boric Acid, 0.4% Sodium chloride, 0.95% Sodium borate

**Purification notes** Ab206195 is purified utilizing a combination of precipitation and/or chromatography techniques.

**Clonality** Polyclonal

**Isotype** IgG

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab206195 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use a concentration of 1 µg/ml.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors