

Product datasheet

Anti-Integrin alpha V beta 3 antibody [27.1 (VNR-1)] ab78289

★★★★★ [5 Abreviews](#) [28 References](#) [1 Image](#)

Overview

Product name	Anti-Integrin alpha V beta 3 antibody [27.1 (VNR-1)]
Description	Mouse monoclonal [27.1 (VNR-1)] to Integrin alpha V beta 3
Host species	Mouse
Tested applications	Suitable for: Flow Cyt Unsuitable for: IHC-P
Species reactivity	Reacts with: Human
Immunogen	The details of the immunogen for this antibody are not available.
General notes	<p>This product was changed from ascites to tissue culture supernatant on 17 May 2019. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	pH: 7.60 Preservative: 0.1% Sodium azide Constituents: 1.45% Sodium chloride, 0.0536% PBS
Purity	Tissue culture supernatant
Purification notes	Purified from TCS.
Clonality	Monoclonal
Clone number	27.1 (VNR-1)

Isotype

IgG1

Applications

The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab78289 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent concentration. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

Application notes

Is unsuitable for IHC-P.

Target


Relevance

Function: The alpha-V (ITGAV) integrins are receptors for vitronectin, cytotactin, fibronectin, fibrinogen, laminin, matrix metalloproteinase-2, osteopontin, osteomodulin, prothrombin, thrombospondin and vWF. They recognize the sequence R-G-D in a wide array of ligands. ITGAV:ITGB3 binds to fractalkine (CX3CL1) and may act as its coreceptor in CX3CR1-dependent fractalkine signaling (PubMed:23125415). ITGAV:ITGB3 binds to NRG1 (via EGF domain) and this binding is essential for NRG1-ERBB signaling (PubMed:20682778). ITGAV:ITGB3 binds to FGF1 and this binding is essential for FGF1 signaling (PubMed:18441324). ITGAV:ITGB3 binds to IGF1 and this binding is essential for IGF1 signaling (PubMed:19578119). ITGAV:ITGB3 binds to PLA2G2A via a site (site 2) which is distinct from the classical ligand-binding site (site 1) and this induces integrin conformational changes and enhanced ligand binding to site 1 (PubMed:18635536, PubMed:25398877). ITGAV:ITGB3 and ITGAV:ITGB6 act as a receptor for fibrillin-1 (FBN1) and mediate R-G-D-dependent cell adhesion to FBN1 (PubMed:12807887, PubMed:17158881). (Microbial infection) Integrin ITGAV:ITGB5 acts as a receptor for adenovirus type C (PubMed:20615244). Integrin ITGAV:ITGB5 and ITGAV:ITGB3 act as receptors for coxsackievirus A9 and B1 (PubMed:9426447, PubMed:15194773, PubMed:7519807). Integrin ITGAV:ITGB3 acts as a receptor for herpes virus 8/HHV-8 (PubMed:18045938). Integrin ITGAV:ITGB6 acts as a receptor for herpes simplex 1/HHV-1 (PubMed:24367260). Integrin ITGAV:ITGB3 acts as a receptor for Human parechovirus 1 (PubMed:11160695). Integrin ITGAV:ITGB3 acts as a receptor for West Nile virus (PubMed:23658209). In case of HIV-1 infection, the interaction with extracellular viral Tat protein seems to enhance angiogenesis in Kaposi's sarcoma lesions (PubMed:10397733). Similarity: Belongs to the integrin alpha chain family. Contains 7 FG-GAP repeats.

Cellular localization

Cell Membrane. Cytoplasmic. Integrin complex. Focal adhesion.

Images



Flow Cytometry - Anti-Integrin alpha V beta 3 antibody [27.1 (VNR-1)] (ab78289)

Overlay histogram showing A431 cells stained with ab78289 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab78289, 2µg/1x10⁶ cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) (**ab96879**) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG1 [ICIGG1] (**ab91353**, 2µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed.

This image was generated using the ascites version of the product.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors