

Product datasheet

Anti-M13 + Fd Bacteriophage Coat Proteins antibody ab6188

★★★★★ [2 Abreviews](#) [10 References](#)

Overview

Product name	Anti-M13 + Fd Bacteriophage Coat Proteins antibody
Description	Rabbit polyclonal to M13 + Fd Bacteriophage Coat Proteins
Host species	Rabbit
Tested applications	Suitable for: Indirect ELISA, ELISA, ICC/IF
Species reactivity	Reacts with: Enterobacteria phage M13
Immunogen	Tissue, cells or virus corresponding to M13 + Fd Bacteriophage Coat Proteins. fd bacteriophage
General notes	Concentration: 5.5 mg/ml- 7.5 mg/ml

Can be used as a capture antibody when coated directly on microtiter plates, or as a primary detection antibody for specifically captured M13 or fd phage. May be useful in rapidly sorting large phage display libraries (antibody, peptide etc.) with the expressed proteins fused to either the gene III or gene VIII protein of the filamentous phage. May be used as a reagent in "phage ELISA" offering sensitive and specific activity for detection of recombinant phages.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	pH: 7.40 Preservative: 0.097% Sodium azide Constituent: 0.0268% PBS
Purity	IgG fraction

Primary antibody notes

Can be used as a capture antibody when coated directly on microtiter plates, or as a primary detection antibody for specifically captured M13 or fd phage. May be useful in rapidly sorting large phage display libraries (antibody, peptide etc.) with the expressed proteins fused to either the gene III or gene VIII protein of the filamentous phage. May be used as a reagent in "phage ELISA" offering sensitive and specific activity for detection of recombinant phages.

Clonality

Polyclonal

Isotype

IgG

Applications**The Abpromise guarantee**

Our **Abpromise guarantee** covers the use of ab6188 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Indirect ELISA		Use at an assay dependent concentration. At least 1/1000 for 5×10^7 phage/ml at least 1/8000 for 5×10^{10} phage/ml.
ELISA	★★★★★ (2)	Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration. PubMed: 17178915

Target**Relevance**

M13 is a filamentous bacteriophage composed of circular single stranded DNA (ssDNA) which is 6407 nucleotides long encapsulated in approximately 2700 copies of the major coat protein P8, and capped with 5 copies of two different minor coat proteins (P9, P6, P3) on the ends.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors