

Anti-Measles Virus Matrix protein antibody ab101015

Overview

Product name	Anti-Measles Virus Matrix protein antibody
Description	Rabbit polyclonal to Measles Virus Matrix protein
Host species	Rabbit
Tested applications	Suitable for: WB, ELISA, ICC
Species reactivity	Reacts with: Measles virus
Immunogen	Synthetic peptides derived from internal parts of Measles Virus Matrix protein.
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Constituent: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab101015 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes

Application	Abreviews	Notes
WB		Use at an assay dependent dilution. Predicted molecular weight: 38 kDa.
ELISA		Use at an assay dependent dilution.
ICC		Use at an assay dependent dilution. (infected cells).

Target

Relevance The Measles Virus Matrix protein has a crucial role in virus assembly and interacts with the RNP complex as well as with the viral membrane.

Cellular localization Envelope protein. Viral matrix protein. Virion

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors