

PE Anti-CD90 / Thy1 antibody [G7] ab24904

★★★★★ [1 Abreviews](#) [8 References](#) [1 Image](#)

Overview

Product name	PE Anti-CD90 / Thy1 antibody [G7]
Description	PE Rat monoclonal [G7] to CD90 / Thy1
Host species	Rat
Conjugation	PE. Ex: 488nm, Em: 575nm
Specificity	ab24904 is specific for Thy 1.1 and Thy 1.2 alloantigens (Thy 1 epitope region A).
Tested applications	Suitable for: Flow Cyt
Species reactivity	Reacts with: Mouse
Immunogen	Tissue, cells or virus corresponding to Mouse CD90/ Thy1. Mouse T cell hybridoma C6/G8
Epitope	Thy 1 epitope region A
General notes	<p>ab24904 stimulates T-cell proliferation and IL2 secretion, via signalling through the T-cell receptor/CD3 complex. ab24904 has also been reported to promote apoptosis of thymocytes and CTL clones and to mediate adhesion of thymocytes to thymic stroma.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Store In the Dark.
Storage buffer	<p>pH: 7.3</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: PBS, 16% Sucrose</p> <p>Stabilising agent.</p>

Purity	Affinity purified
Primary antibody notes	ab24904 stimulates T-cell proliferation and IL2 secretion, via signalling through the T-cell receptor/CD3 complex. ab24904 has also been reported to promote apoptosis of thymocytes and CTL clones and to mediate adhesion of thymocytes to thymic stroma.
Clonality	Monoclonal
Clone number	G7
Isotype	IgG2c
Light chain type	kappa

Applications


The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab24904 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt	★★★★★ (1)	Use 2µg for 10 ⁶ cells. ab154450 - Rat monoclonal IgG2c, is suitable for use as an isotype control with this antibody.

Target

Function	May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.
Sequence similarities	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
Cellular localization	Cell membrane.

Images



Flow cytometry analysis staining CD90 / Thy1 in mouse splenocytes using ab24904 at a dilution of 0.03 µg/10⁶ cells.

Flow Cytometry - PE Anti-CD90 / Thy1 antibody
[G7] (ab24904)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors