

Product datasheet

PerCP Anti-CD8 antibody [MEM-31] ab65949

★★★★★ [1 Abreviews](#) [2 References](#)

Overview

Product name	PerCP Anti-CD8 antibody [MEM-31]
Description	PerCP Mouse monoclonal [MEM-31] to CD8
Host species	Mouse
Conjugation	PerCP. Ex: 482nm, Em: 675nm
Specificity	ab65949 recognizes a conformationally-dependent epitope of CD8. The antibody does not react with formaldehyde-fixed cells and is negative in Western Blotting.
Tested applications	Suitable for: Flow Cyt Unsuitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Tissue, cells or virus corresponding to Human CD8. Crude thymus membrane fraction
Positive control	human blood
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.097% Sodium azide Constituents: PBS, 0.2% BSA
Purity	Size exclusion
Purification notes	>95% by SDS-PAGE
Clonality	Monoclonal
Clone number	MEM-31
Isotype	IgG2a

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab65949 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt	★★★★★ (1)	Use 10µl for 10 ⁶ cells. ab110422 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody. Use Neat: 10µl for 10 ⁶ cells, or 100µl whole blood. Is unsuitable for applications using formaldehyde fixed cells.

Application notes Is unsuitable for WB.

Target

Function	Identifies cytotoxic/suppressor T-cells that interact with MHC class I bearing targets. CD8 is thought to play a role in the process of T-cell mediated killing. CD8 alpha chains binds to class I MHC molecules alpha-3 domains.
Involvement in disease	Defects in CD8A are a cause of familial CD8 deficiency (CD8 deficiency) [MIM:608957]. Familial CD8 deficiency is a novel autosomal recessive immunologic defect characterized by absence of CD8+ cells, leading to recurrent bacterial infections.
Sequence similarities	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
Post-translational modifications	All of the five most C-terminal cysteines form inter-chain disulfide bonds in dimers and higher multimers, while the four N-terminal cysteines do not.
Cellular localization	Secreted and Cell membrane.
Form	CD8 beta tissue specificity: Isoform 1, isoform 3, isoform 5, isoform 6, isoform 7 and isoform 8 are expressed in both thymus and peripheral CD8+ T-cells. Expression of isoform 1 is higher in thymus CD8+ T-cells than in peripheral CD8+ T-cells. Expression of isoform 6 is higher in peripheral CD8+ T-cells than in thymus CD8+ T-cells. CD8 beta PTM: Phosphorylated as a consequence of T-cell activation.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors