

Product datasheet

Anti-Respiratory Syncytial Virus Nucleoprotein antibody [B325M] ab252771

Overview

Product name	Anti-Respiratory Syncytial Virus Nucleoprotein antibody [B325M]
Description	Mouse monoclonal [B325M] to Respiratory Syncytial Virus Nucleoprotein
Host species	Mouse
Specificity	RSV Nucleoprotein (42 kDa).
Immunogen	Tissue, cells or virus corresponding to Respiratory Syncytial Virus Nucleoprotein. RSV strains 127, SNK and 9007, Long, Randall, 8/60 and A2.
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.0975% Sodium azide Constituent: PBS
Purity	Protein A purified
Purification notes	>90% pure (SDS-PAGE). Sourced from tissue culture supernatant.
Clonality	Monoclonal
Clone number	B325M
Isotype	IgG1

Target

Relevance	Respiratory syncytial virus (RSV) is a major cause of respiratory illness in young children. RSV infection produces a variety of signs and symptoms involving different areas of the respiratory tract, from the nose to the lungs. RSV is a negative sense, enveloped RNA virus. The virion is variable in shape and size with average diameter of between 120 and 300 nm.
Cellular localization	Virus particle

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors