

Product datasheet

Anti-Sox9 α antibody - N-terminal ab209820

[1 References](#) [1 Image](#)

Overview

Product name	Anti-Sox9a antibody - N-terminal
Description	Rabbit polyclonal to Sox9a - N-terminal
Host species	Rabbit
Tested applications	Suitable for: IHC - Wholemout
Species reactivity	Reacts with: Zebrafish
Immunogen	Recombinant fragment within Zebrafish Sox9a aa 20-201 (N terminal). The exact sequence is proprietary. Database link: Q9DFH2
Positive control	Zebrafish embryos.
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 79.99% PBS, 20% Glycerol (glycerin, glycerine)
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab209820 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC - Wholemout		1/100 - 1/500.

Target

Images



Immunohistochemical analysis (whole mount) of paraformaldehyde-fixed Zebrafish embryo labeling Sox9a with ab209820 at 1/200 dilution.

IHC - Wholemout - Anti-Sox9a antibody - N-terminal (ab209820)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors