

Texas Red ® Anti-Biotin antibody ab6653

1 Image

Overview

Product name	Texas Red ® Anti-Biotin antibody
Description	Texas Red ® Goat polyclonal to Biotin
Host species	Goat
Conjugation	Texas Red ®. Ex: 596nm, Em: 620nm
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Species independent
Immunogen	Chemical/ Small Molecule conjugated to KLH.
General notes	<p>Texas Red™ is a registered trademark of Molecular Probes Inc.</p> <p>Texas Red™ Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 2.5 moles Texas Red™ per mole of Goat IgG.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Biotin coupled to sepharose beads.

Clonality	Polyclonal
Isotype	IgG

Applications


The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab6653 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration.

Target

Relevance Biotin is a water soluble vitamin, generally classified as a B complex vitamin, also called vitamin B4. After the initial discovery of biotin, nearly forty years of research were required to establish it as a vitamin. Biotin is required by all organisms but can only be synthesized by bacteria, yeasts, molds, algae, and some plant species. Biotin is required as prosthetic group of enzymes involved in incorporation of carbon dioxide into organic compounds. Biotin has a MW of 244 Da.

Images



Texas Red ® Anti-Biotin antibody (ab6653) at 1/1000 dilution (60 minutes at room temperature) + Biotin conjugated Guinea Pig IgG separated by SDS-PAGE at 0.05 µg

Observed band size: 25,55 kDa

Blocking Buffer for Fluorescent Western Blotting for 30 min at room temperature.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors