

Product datasheet

Anti-Yellow Fever Virus antibody [3576] ab22839

★★★★★ 1 Abreviews 1 References

Overview

Product name	Anti-Yellow Fever Virus antibody [3576]
Description	Mouse monoclonal [3576] to Yellow Fever Virus
Host species	Mouse
Specificity	Does not detect other Flaviviruses.
Tested applications	Suitable for: ELISA
Species reactivity	Reacts with: Yellow fever virus
Immunogen	Tissue, cells or virus corresponding to Yellow fever virus Yellow Fever Virus. Whole Yellow fever virus preparation.
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	<p>pH: 7.20</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituent: 0.0268% PBS</p>
Purity	Protein A purified
Clonality	Monoclonal
Clone number	3576
Isotype	IgG2a

Applications

The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab22839 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/20 - 1/200.

Target

Relevance

Yellow Fever Virus is a mosquito-borne species of the genus Flavivirus that causes yellow fever in Central and South America and Africa. Yellow fever (also called black vomit or sometimes The American Plague) is an acute viral disease. It is still an important cause of haemorrhage illness in several African and South American countries despite existence of an effective vaccine. In the past it was a source of several devastating epidemics.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors