

Product datasheet

Anti-ZNF545 antibody [OTI5D8] ab279394

1 Image

Overview

Product name	Anti-ZNF545 antibody [OTI5D8]
Description	Mouse monoclonal [OTI5D8] to ZNF545
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Recombinant full length protein corresponding to Human ZNF545 aa 1-532. NP_597723 Database link: Q8N141
Positive control	WB: HepG2 cell extracts.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.3 Preservative: 0.02% Sodium azide Constituents: 1% BSA, 50% Glycerol (glycerin, glycerine)
Purity	Protein A/G purified
Purification notes	Purified from TCS.
Clonality	Monoclonal
Clone number	OTI5D8
Isotype	IgG2a

Applications


The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab279394 in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/2000. Predicted molecular weight: 63 kDa.

Target

Function	May be involved in transcriptional regulation.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 13 C2H2-type zinc fingers. Contains 1 KRAB domain.
Cellular localization	Nucleus.

Images



Anti-ZNF545 antibody [OTI5D8] (ab279394) at 1/500 dilution +
HepG2 (Human liver hepatocellular carcinoma cell line) cell extracts
at 35 µg

Predicted band size: 63 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors