

Human Lano peptide ab39803

Description

Product name	Human Lano peptide
Animal free	No
Nature	Synthetic
Species	Human

Specifications

Our **Abpromise guarantee** covers the use of **ab39803** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot Blocking
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.20 Preservative: 0.02% Sodium azide Constituents: PBS, 0.1% BSA
------------------------------	--

General Info

Tissue specificity	Expressed strongly in testis and placenta, followed by heart, lung, kidney, thyroid, trachea, colon, prostate and pancreas.
Sequence similarities	Contains 17 LRR (leucine-rich) repeats.
Cellular localization	Cytoplasm. Membrane. Localized at the basolateral side of epithelial cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors