

Product datasheet

Native Human Chylomicrons protein ab285812

[1 Image](#)

Description

Product name	Native Human Chylomicrons protein
Purity	<p>>= 95 % SDS-PAGE.</p> <p>Single arc by IEP against antisera to whole human serum. Essentially free of other plasma lipoproteins as determined by electrophoresis using a SPIFE Vis Cholesterol gel kit for lipids and Coomassie Blue for proteins. Content is 98% lipid, 2% protein. Density is <0.95 g/ml</p>
Protein length	Full length protein
Animal free	No
Nature	Native
Species	Human
Molecular weight information	50-1000 kDa molecular weight
Additional sequence information	Fresh, non-frozen Human Plasma. Prepared from plasma shown to be non-reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests.

Specifications

Our **Abpromise guarantee** covers the use of **ab285812** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid
Additional notes	<p>This product is manufactured by BioVision, an Abcam company and was previously called 7285 Chylomicrons, Human Plasma. 7285-1000 is the same size as the 1 mg size of ab285812.</p> <p>Chylomicrons are the large triglyceride-rich lipoproteins of enteric origin that carry exogenous triglycerides into the plasma. Circulation of chylomicrons distributes dietary triglyceride to adipose tissue and to muscle. Found in post prandial plasma, chylomicrons generally clear the blood stream in a few hours.</p>

Preparation and Storage

Stability and Storage	<p>Shipped at 4°C. Store at +4°C. Do Not Freeze. Store In the Dark.</p> <p>pH: 7.40</p> <p>Constituents: 0.877% Sodium chloride, 0.01% EDTA</p>
------------------------------	---

Images



12% BIO-RAD Tris HCl gel

Lane 1. MW Standard

Lane 2. Chylomicrons - 10 µg (reduced / heat)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors