

Product datasheet

Recombinant Dog Interferon gamma protein ab222448

1 Image

Description

Product name	Recombinant Dog Interferon gamma protein	
Purity	> 95 % SDS-PAGE.	
Endotoxin level	< 1.000 Eu/μg	
Expression system	Escherichia coli	
Accession	<u>P42161</u>	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Dog	
Sequence	MQAMFFKEIENLKEYFNASNPDVSDGGSLFVDILKKWRE ESDKTIQSQI VSFYKLFDFNFKDNQIQRSMDTIKEDMLGKFLNSSTSKRE DFLKLIQIP VNDLQVQRKAINELIKVMNDLSPRSNLRKRKRSQNLFRGR RASK	
Predicted molecular weight	17 kDa	
Amino acids	24 to 166	
Additional sequence information	This product is the mature full length protein from aa 24 to 166. The signal peptide is not included.	

Specifications

Our **Abpromise guarantee** covers the use of **ab222448** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 10% Glycerol (glycerin, glycerine), PBS

General Info

Function

Produced by lymphocytes activated by specific antigens or mitogens. IFN-gamma, in addition to having antiviral activity, has important immunoregulatory functions. It is a potent activator of macrophages, it has antiproliferative effects on transformed cells and it can potentiate the antiviral and antitumor effects of the type I interferons.

Tissue specificity

Released primarily from activated T lymphocytes.

Involvement in disease

In Caucasians, genetic variation in IFNG is associated with the risk of aplastic anemia (AA) [MIM:609135]. AA is a rare disease in which the reduction of the circulating blood cells results from damage to the stem cell pool in bone marrow. In most patients, the stem cell lesion is caused by an autoimmune attack. T-lymphocytes, activated by an endogenous or exogenous, and most often unknown antigenic stimulus, secrete cytokines, including IFN-gamma, which would in turn be able to suppress hematopoiesis.

Sequence similarities

Belongs to the type II (or gamma) interferon family.


Post-translational modifications

Proteolytic processing produces C-terminal heterogeneity, with proteins ending alternatively at Gly-150, Met-157 or Gly-161.

Cellular localization

Secreted.

Images



3 µg ab222448 analyzed by 15% SDS-PAGE.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors