

Product datasheet

Recombinant Fimbrial subunit type 3 protein (Tagged) ab236181

1 Image

Description

| | |
|--|---|
| Product name | Recombinant Fimbrial subunit type 3 protein (Tagged) |
| Purity | > 90 % SDS-PAGE. |
| Expression system | Escherichia coli |
| Accession | <u>P12267</u> |
| Protein length | Full length protein |
| Animal free | No |
| Nature | Recombinant |
| Sequence | ADTNVGGGQVNFFGKVTDVSCTVSVNGQGSDANVYLSP VTLTEVKAAAAD TYLKPKSFTIDVSDCQAADGTKQDDVSKLGVNWTGGNLL AGATAKQQGYL ANTEAAGAQNIQLVLSTDNATALTNKIIPGDSTQPKAAGDA SAVQDGARF TYYVGYATSTPTTVTTGGVNSYATYEITYQ |
| Predicted molecular weight | 35 kDa including tags |
| Amino acids | 23 to 202 |
| Tags | His tag N-Terminus |
| Additional sequence information | Klebsiella pneumoniae. This product is the mature full length protein from aa 23 to 202. The signal peptide is not included. N-terminal 6xHis-SUMO-tag. |

Specifications

Our **Abpromise guarantee** covers the use of **ab236181** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid


Preparation and Storage

Stability and Storage Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.2

Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) analysis with 5% enrichment gel and 15% separation gel of ab236181.

SDS-PAGE - Recombinant Fimbrial subunit type 3 protein (Tagged) (ab236181)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors