

Product datasheet

Recombinant Hemagglutinin neuraminidase protein (Tagged) ab236804

1 Image

Description

Product name	Recombinant Hemagglutinin neuraminidase protein (Tagged)
Purity	> 85 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>P35741</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Sequence	NGAANNSGCGAPVHDPDYGGIGKELVDDASDVTSFYPS AFQEHLNFIP APTTGSGCTRIPSFDISATHYCYTHNVILSGCRDHSHSHQY LALGVLRTS ATGRVFFSTLRSINLDDNQRKSCSVSATPLGCDMLCSKI TETEEEDYSS VTPTSMVHGRLGFDGQYHEKDLDVITLFDKDWVANYPGVG GGSFIDNRVWF PVYGGGLKPNSPSTVQEGRYVIYKRYNDTCPDEQDYQIRM AKSSYKPGRF GGKRVQQAILSIVSTSLGEDPVLTIPTNTVTLMGAEGRVL TVGTSHFLY QRGSSYFSPALLYPMTVNNKTATLHSPYTFNAFTRPGSVP CQASARCPNS CVTGVYTD
Predicted molecular weight	44 kDa including tags
Amino acids	115 to 473
Tags	His tag N-Terminus
Additional sequence information	Newcastle disease virus (strain Her/33). N-terminal 10xHis-tagged and C-terminal Myc-tagged.

Specifications

Our **Abpromise guarantee** covers the use of **ab236804** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage


Stability and Storage

Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.20

Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)

Images



ab236804 analyzed by (Tris-Glycine gel) discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors