

Product datasheet

Recombinant HIV-1 p31 protein ab173265

Description

Product name	Recombinant HIV-1 p31 protein
Purity	> 95 % SDS-PAGE. ab173265 is purified using S-Sepharose, Ceramic Hydroxyapatite and S-300.
Expression system	Escherichia coli
Protein length	Full length protein
Animal free	No
Nature	Recombinant

Specifications

Our **Abpromise guarantee** covers the use of **ab173265** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot
Form	Liquid
Additional notes	ab173265 is immunoreactive with all sera of HIV-infected individuals with immuno response to the HIV POL. The protein is a GST N-terminal fusion.

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. pH: 8.00 Constituents: 50% Glycerol, 0.61% Tris
------------------------------	---

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors