

Product datasheet

Recombinant Human ACSM3/SAH protein ab127250

Description

Product name	Recombinant Human ACSM3/SAH protein
Purity	> 80 % SDS-PAGE. Purified via His tag
Expression system	Escherichia coli
Accession	<u>Q53FZ2</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Predicted molecular weight	26 kDa
Amino acids	36 to 264
Tags	His-DHFR tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab127250** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Lyophilized
Additional notes	This product was previously labelled as ACSM3

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
Reconstitution	Reconstitute with water to desired concentration.

General Info

Function	Has medium-chain fatty acid:CoA ligase activity with broad substrate specificity (in vitro). Acts on acids from C(4) to C(11) and on the corresponding 3-hydroxy- and 2,3- or 3,4-unsaturated acids (in vitro).
Sequence similarities	Belongs to the ATP-dependent AMP-binding enzyme family.
Cellular localization	Mitochondrion matrix.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors