

Recombinant Human ADH1B protein ab116934

1 References 1 Image

Description

Product name	Recombinant Human ADH1B protein		
Expression system	Wheat germ		
Accession	P00325		
Protein length	Full length protein		
Animal free	No		
Nature	Recombinant		
Species	Human		
Sequence	MSTAGKVIKCKAAVLWEVKKPFSIEDVEVAPPKAYEVRIK MVAVGICRTD DHVVSGNLVTPLPVILGHEAAGVESVGEGVTTVKPGDKVI PLFTPQCGK CRVCKNPESNYCLKNDLG NPRGTLQDGTRRFTCRGKPIH HFLGTSTFSQY TVVDENAVAKIDAASPLEKVCLIGCGFSTGYGSAVNVAKV TPGSTCAVFG LGGVGLSAVMGCKAAGAARIAVDINKDKFAKAKELGATE CINPQDYKKP IQEVLKEMTDGGVDFSFEVIGRLDTMMASLLCCHEACGTS VIVGVPPASQ NLSINPMLLLTGRTWKGAVYGGFKSKEGIPKLVADFMAKK FSLDALITHV LPFEKINEGFDLLHSGKSIRTVLTF		
Predicted molecular weight	67 kDa including tags		
Amino acids	1 to 375		

Specifications

Our Abpromise guarantee covers the use of ab116934 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	ELISA SDS-PAGE Western blot
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

General Info


Sequence similarities

Belongs to the zinc-containing alcohol dehydrogenase family.

Cellular localization

Cytoplasm.

Images



12.5% SDS-PAGE stained with Coomassie Blue showing ab116934 at approximately 66.99 kDa.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors