

Recombinant Human AFM protein ab114886

1 Image

Description

Product name	Recombinant Human AFM protein
Expression system	Wheat germ
Accession	<u>P43652</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MKLLKLTGFIFFLFFLTESLTLPTQPRDIENFNSTQKFIEDNI EYITIAFAQYVQEATF EEMEKLVKDMVEYKDRCMADKTLPECSKLPNNVLQEKI CAMEGLPQKHNFSHCCSKVDAQ RRLCFFYNKKSDVGFLPPFPPTLDPEE KCQAYESNRESLLNHFLYEVARRNPFVFAPTLT VAVHFEEVAKSCCE EQNKVNCLQTRAIPVTQYLKAFSSYQKHVCGALLKFGTKV VHFYAI LSQKFPKIEFKELISLVEDVSSNYDGCCEGDVVQCIRDTS KVMNHICKSQDSISSKIK ECCEKKIPERGQCINSNKDDRPKDLSLREGKFTDSENV QERDADPDFTFAKFTFEYSR RHPDLSIPELLRIVQYKDLLRNCCNTE NPPGCYRYAEDKFNETTEKSLKMOVQEQECKHFQ NLGKDGLKYHYLIRLT KIAPQLSTEELVSLGEKMTAFTTCCTLSEEFACVDNLAD LVFGEELCG VNENRTINPAVDHCCKTNFAFRRPCFESLKADKTYPPPF SQDLFTFHADMCSQN EELQRKTDRFLVNLVKLKHETDEELQSLFTNFANVVDKC CKAESPEVCFNEESPKIGN
Predicted molecular weight	96 kDa including tags
Amino acids	1 to 599

Specifications

Our **Abpromise guarantee** covers the use of **ab114886** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	SDS-PAGE
	Western blot
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
	pH: 8.00
	Constituents: 0.3% Glutathione, 0.79% Tris HCl

General Info

Function	Vitamin E binding protein. May transport vitamin E in body fluids under conditions where the lipoprotein system is not sufficient or across the blood-brain barrier.
Sequence similarities	Belongs to the ALB/AFP/VDB family. Contains 3 albumin domains.
Post-translational modifications	N-glycosylated.
Cellular localization	Secreted.

Images



12.5% SDS-PAGE showing ab114886 at approximately 96.10kDa.
Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors