

Product datasheet

Recombinant Human AHNAK protein ab163959

[1 Image](#)

Description

Product name	Recombinant Human AHNAK protein
Expression system	Wheat germ
Accession	<u>Q09666-2</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MEKEETTRELLLPNWQSGSGSHGLTIAQRDDGVFVQEVTQ NSPAARTGVVK EGDQIVGATYFDNLQSGEVTQLLNTMGHHTVGLKLHRKG DRSPEPGQTW TREVFSSCSSEVVLNTPQPSALECKDQNKQKEASSQAG AVSVSTPNAGL
Amino acids	1 to 149
Tags	GST tag N-Terminus
Additional sequence information	Isoform 2 of AHNAK.

Specifications

Our **Abpromise guarantee** covers the use of **ab163959** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot

Form	Liquid
-------------	--------

Additional notes

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------------------------	---


pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function	May be required for neuronal cell differentiation.
Sequence similarities	Contains 1 PDZ (DHR) domain.
Cellular localization	Nucleus.

Images



SDS-PAGE - Recombinant Human AHNAK protein
(ab163959)

ab163959 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors