

Recombinant Human APOBEC4 protein (denatured)
ab171472

1 Image

Description

Product name	Recombinant Human APOBEC4 protein (denatured)
Purity	> 85 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>Q8WW27</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMEPIYEE YLANHGTVK PYYWLSFSLD CSNCPYHIRT GEEARVSLTE FCQIFGFPYG TTFPQTKHLT FYELKTSSGS LVQKGHASSC TGNYHPESM LFEMNGYLDs AIYNDSIRH IILYSNNSPC NEANHCCISK MYNFLITYPG ILSIYFSQL YHTEMDFPAS AWNREALRSL ASLWPRVVLS PISGGWHSV LHSFISGVSG SHVFQPILTG RALADRHNA Y EINATGVKP YFTDVLLQTK RNPNTKAQEA LESYPLNNAF PGQFFQMPSG QLQPNLPPDL RAPVVFVLVP LRDLPPMHMG QNPKNPRNIV RHLNMPQMSF QETKDLGRLP TGRSVEVEI TEQFASSKEA DEKKKKKGKK
Predicted molecular weight	44 kDa including tags
Amino acids	1 to 367
Tags	His tag N-Terminus
Description	Recombinant Human APOBEC4 protein

Specifications

Our **Abpromise guarantee** covers the use of **ab171472** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 2.4% Urea, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine)


General Info

Function Putative C to U editing enzyme whose physiological substrate is not yet known.

Tissue specificity Predominantly expressed in testis.

Sequence similarities Belongs to the cytidine and deoxycytidylate deaminase family.

Images



15% SDS-PAGE analysis of 3µg ab171472.

SDS-PAGE - Recombinant Human APOBEC4
protein (denatured) (ab171472)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors