

Recombinant Human ASC1 protein ab152970

1 Image

Description

Product name	Recombinant Human ASC1 protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MAVAGAVSGEPLVHWCTQQLRKTFGLDVSEEIIQYVLSIE SAEEIREYVT DLLQGNEGKKKGQFIEELITKWQKNDQELISDPLQQCFKKD EILDGQKSGD HLKRGRKKGRNRQEVPAFTEPDTTAEVKTPFDLAKAQEN SNSVKKKTKFV NLYTREGQDRLAVLLPGRHPCDCLGQKHKLINNCLICGRIV CEQEGSGPC LFCGTLVCTHEEQDILQRDSNKSQKLLKKLMSGVENS GK VDISTKDLLPH QELRIKSGLEKAIKHKDKLLEFDRTSIRRTQVIDDESDYFAS DSNQWLSK LERETLQKREEELRELRHASRLSKKVTIDFAGRKILEEENS LAEYHSRLD ETIQAIANGTLNQPLTKLDRSSEEPLGVLVNP NMYQSPPQ WVDHTGAASQ KKA FRSSGFGLEFNSFQHQLRIQDQEFQEGFDGGWCLS VHQPWASLLVRG IKRVEGRSWYTPHRGRLWIAATAKKPSPQEVSELQATYRL LRGKDVEFPN DYPSGCLLGCVDLIDCLSQKQFKEQFPDISQESDSPFVFI CKNPQEMVVK FPIKGNPKIMWKLD SKIHQGAKKGLMKQNKAV
Amino acids	1 to 581
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab152970** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Western blot

ELISA

Form Liquid

Additional notes

Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00


Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function Transcription coactivator of nuclear receptors which functions in conjunction with CBP-p300 and SRC-1 and may play an important role in establishing distinct coactivator complexes under different cellular conditions. Plays a pivotal role in the transactivation of NF-kappa-B, SRF and AP1. Acts as a mediator of transrepression between nuclear receptor and either AP1 or NF-kappa-B. Plays a role in androgen receptor transactivation and in testicular function.

Cellular localization Nucleus. Cytoplasm. Cytoplasmic under conditions of serum deprivation.

Images



ab152970 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors