

Recombinant Human BCL6B protein ab166053

1 Image

Description

Product name	Recombinant Human BCL6B protein		
Expression system	Wheat germ		
Accession	<u>Q8N143</u>		
Protein length	Full length protein		
Animal free	No		
Nature	Recombinant		
Species	Human		
Sequence	MGSPAAPEGALGYVREFTRHSSDVLGNLNELRLRGILTDV TLLVGGQPLR AHKAVLIACSGFFYSIFRGRAGVGVDVLSLPGGPEARGFA PLLDIFYTSR LRLSPATAPAVLAAATYLQMEHVQACHRFIQASYEPLGIS LRPLEAEP TPPTAPPPGSPRRSEGHDPPTESRSCSQGPPSPASPD PKACNWKKYKI VLNSQASQAGSLVGERSSGQPCQARLPSGDEASSSS SSSSSSSEEGPI PGPQSRLSPTAATVQFKCGAPASTPYLLTSQAQDTSGSP SERARPLPGSE FFSCQNSEAVAGCSSGLDSLVPGEDDKPYKCQLCRSSF RYKGNLASHRTV HTGEKPYHCSICGARFNRPANLKTHSRIHSGEKPYKCETC GSRFVQVAHL RAHVLIHTGEKPYPCPTCGTRFRHLQTLKSHVRIHTGEKPY HCDPCGLHF RHKSQLRLHLRQKHGAATNTKVHYHILGGP		
Amino acids	1 to 480		
Tags	GST tag N-Terminus		

Specifications

Our **Abpromise guarantee** covers the use of **ab166053** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot
--------------	--------------

ELISA

Form Liquid

Additional notes

Preparation and Storage


Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Relevance BCL6B acts as a sequence-specific transcriptional repressor in association with BCL6. It may function in a narrow stage or be related to some events in the early B-cell development.

Cellular localization Nuclear

Images



SDS-PAGE - Recombinant Human BCL6B protein
(ab166053)

ab166053 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors