

## Product datasheet

# Recombinant Human BFAR protein ab162491

1 Image

### Description

<b>Product name</b>	Recombinant Human BFAR protein
<b>Expression system</b>	Wheat germ
<b>Protein length</b>	Full length protein
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Sequence</b>	<pre> MEEPQKSYVNTMDLERDEPLKSTGPGQISVSEFSCHCCYDI LVNPTTLNCG HSFCRHCLALWWASSKKTECPECREKWEGFPKVSILLR DAIEKLFDAIR LRFEDIQQNNDIVQSLAAFQKYGNDQIPLAPNTGRANQQM GGGFFSGVLT ALTGVAVVLLVYHWSSRESEHDLLVHKAVAKWTAEVVL WLEQLGPWASL YRERFLSERVNGRLLLTLTEEFKSKTPYTIENSSHRRAILME LERVKALG VKPPQNLWEYKAVNPGRSFLLYALKSSPRLSLLYLFD YDTFLPFIH TICPLQEDSSGEDIVTKLLDLKEPTWKQWREFLVKYSFLP YQLIAEFAWD WLEVHYWTSRFLINAMLLSVLELFSFWRWSRSELKTPVQ RMWSHFVKV STQGLFVAMFWPLIPQFVCNCLFYWALYFNPIINDLVVKEL RRLETQVL </pre>
<b>Amino acids</b>	1 to 450
<b>Tags</b>	GST tag N-Terminus

### Specifications

Our **Abpromise guarantee** covers the use of **ab162491** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	Western blot
	ELISA

**Form** Liquid

**Additional notes**

---

**Preparation and Storage**

**Stability and Storage** Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.  
pH: 8.00  
Constituents: 0.31% Glutathione, 0.79% Tris HCl

---

**General Info**

**Function** Apoptosis regulator. Has anti-apoptotic activity, both for apoptosis triggered via death-receptors and via mitochondrial factors.

**Tissue specificity** Expressed highly in brain, moderately in small intestine, weakly in testes and only faintly in liver and skeletal muscle. Not expressed in heart, kidney, lung and spleen.

**Sequence similarities** Contains 1 RING-type zinc finger.  
Contains 1 SAM (sterile alpha motif) domain.

**Cellular localization** Endoplasmic reticulum membrane.

---

**Images**

---



ab162491 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human BFAR protein  
(ab162491)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

---

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors