

Product datasheet

Recombinant Human Cathelicidin/CLP protein (denatured) **ab140725**

[2 References](#) [1 Image](#)

Description

Product name	Recombinant Human Cathelicidin/CLP protein (denatured)
Purity	> 85 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>P49913</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSQVLSYKE AVLRAIDGIN QRSSDANLYR LLDLDPRPTM DGDPTDPKPV SFTVKETVCP RTTQQSPEDC DFKKDGLVKR CMGTVTLNQA RGSFDISCDK DNKRFALLGD FFRKSKEKIG KEFKRIVQRI KDFLRNLVPR TES
Predicted molecular weight	18 kDa including tags
Amino acids	30 to 170
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab140725** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid
Additional notes	This product was previously labelled as Cathelicidin

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 2.4% Urea, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine)

General Info

Function	Binds to bacterial lipopolysaccharides (LPS), has antibacterial activity.
Tissue specificity	Expressed in bone marrow and testis and neutrophils.
Sequence similarities	Belongs to the cathelicidin family.
Post-translational modifications	The N-terminus is blocked.
Cellular localization	Secreted.

Images



SDS-PAGE analysis of ab140725 (3 µg) under reducing conditions and visualized by coomassie blue stain.

SDS-PAGE - Recombinant Human Cathelicidin/CLP protein (denatured) (ab140725)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors