

Product datasheet

Recombinant Human CBF1 interacting corepressor/CIR1 protein ab160540

1 Image

Description

Product name	Recombinant Human CBF1 interacting corepressor/CIR1 protein		
Expression system	Wheat germ		
Accession	<u>Q86X95</u>		
Protein length	Full length protein		
Animal free	No		
Nature	Recombinant		
Species	Human		
Sequence	MGKSFANFMCKKDFHPASKSNIKKVWMAEQKISYDKKKQ EELMQQYLKEQ ESYDNRLLMGDERVKNGLNFMYEAPPGAKKENKEKEETE GETEYKFEWQK GAPREKYAKDDMNIRDQPFGIQVRNVRCIKCHKWGHVNT DRECPLFGLSG INASSVPTDGSGPSMHPSELIAEMRNSGFALKRNVLGRNL TANDPSQEYV ASEGEEDPEVEFLKSLTTKQKQKLLRKLDRLKKKKKKKD RKKKKKFQKSRS KHKKHKSSSSSSSSSSSSSTETSESSSESESNNKEKKI QRKKRKKNKCS GHNNSDSEEKDKSKKRKLHEELSSSHHNREKAKEKPRFL KHESREDSKW SHSDSDKKSRTHKHSPEKRGSERKEGSSRSHGREERSR RSRSRSPGSYKQ RETRKRAQRNPGEEQSRRNDSRSHGTDLYRGEKMYREH PGGTHTKVTQRE		
Amino acids	1 to 450		
Tags	GST tag N-Terminus		

Specifications

Our **Abpromise guarantee** covers the use of **ab160540** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot
Form	Liquid
Additional notes	This product was previously labelled as CBF1 interacting corepressor.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
	pH: 8.00
	Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Relevance	CBF1 interacting corepressor may modulate splice site selection during alternative splicing of pre-mRNAs. It regulates transcription and acts as corepressor for RBPSUH; recruits RBPSUH to the Sin3-histone deacetylase complex (HDAC).
Cellular localization	Nuclear

Images



ab160540 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors