

### Recombinant Human CD62L protein ab200299

#### Description

<b>Product name</b>	Recombinant Human CD62L protein
<b>Purity</b>	> 90 % SDS-PAGE.
<b>Expression system</b>	Escherichia coli
<b>Accession</b>	<b><u>P14151</u></b>
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Sequence</b>	<p>MASMTGGQQMGRGHHHHHHGNLYFQGEFLCCDFLAHHG  TDCWTYHYSEKP  MNWQRARRFCRDNYDLVAIQNKAEIEYLEKTLPFSSRSYY  WIGIRKIGGI  WTWVGTKNSL TEEAENWGDGEPNNKKNKEDCVEIMKRN  KDAGKWND DAC  HKLKAALCYTASCQPWSCSGHGECVEIINNYTCNCDVGYY  GPQCQFVIQC  EPLAEPGLTMDCTHPLGNFSFSSQCAFSCSEGTNLTGIE  ETTCGPFGNW  SSPEPTCQVIQCEPLSAPDLGIMNC SHPLASFSFTSACTFI  CSEGT ELIG  KKKTICESSGWSNPSPICQKLDKSF SMIKEGDYN</p>
<b>Predicted molecular weight</b>	38 kDa including tags
<b>Amino acids</b>	29 to 332
<b>Tags</b>	His-T7 tag N-Terminus
<b>Additional sequence information</b>	Human CD62L extracellular domain

#### Specifications

Our **Abpromise guarantee** covers the use of **ab200299** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** SDS-PAGE

**Form** Liquid

## Preparation and Storage

---

**Stability and Storage** Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituent: 0.32% Tris HCl

Contains NaCl, KCl, EDTA, arginine, DTT and Glycerol.

## General Info

---

**Function** Cell surface adhesion protein. Mediates the adherence of lymphocytes to endothelial cells of high endothelial venules in peripheral lymph nodes. Promotes initial tethering and rolling of leukocytes in endothelia.

**Tissue specificity** Expressed in B cell lines and T lymphocytes.

**Sequence similarities** Belongs to the selectin/LECAM family.  
Contains 1 C-type lectin domain.  
Contains 1 EGF-like domain.  
Contains 2 Sushi (CCP/SCR) domains.

**Cellular localization** Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors