

Product datasheet

Recombinant Human Chemerin protein ab256228

1 Image

Description

Product name	Recombinant Human Chemerin protein
Purity	> 95 % SDS-PAGE.
Endotoxin level	< 1.000 Eu/μg
Expression system	Escherichia coli
Accession	<u>Q99969</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MELTEAQRRLQVALEEFHKHPPVQWAFQETSVESA VDT TPFPAGIFVRLE FKLQQTSCRKRDWKKPECKVRPNRKRKCLACIKLGSE DKVLGRLVHCPI ETQVLREAEEHQETQCLRVQRAGEDPHSFYFPGQFAFS
Predicted molecular weight	16 kDa
Amino acids	21 to 157
Additional sequence information	Full length chain without signal peptide, without propeptide.

Specifications

Our **Abpromise guarantee** covers the use of **ab256228** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Lyophilized

Preparation and Storage

Stability and Storage	Shipped at room temperature. Store at -20°C. Constituent: 0.1% Trifluoroacetic acid
------------------------------	--

0.2 micron filtered

Reconstitution

Sterile water to 0.1 mg/mL. Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.1% BSA solution, store at -80°C and avoid repeat freeze thaws.

General Info


Tissue specificity

Highly expressed in skin (basal and suprabasal layers of the epidermis, hair follicles and endothelial cells). Also found in pancreas, liver, spleen, prostate, ovary, small intestine and colon.

Cellular localization

Secreted.

Images



4-20% SDS-PAGE - Recombinant Human Chemerin protein (ab256228 1µg per lane).

Lane 1: Non-reducing conditions.

Lane 2: Reducing conditions.

Coomassie blue staining.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors