

## Product datasheet

# Recombinant Human Chromogranin B protein ab93930

[1 Image](#)

### Description

<b>Product name</b>	Recombinant Human Chromogranin B protein
<b>Purity</b>	> 80 % SDS-PAGE. Purified using conventional chromatography techniques
<b>Expression system</b>	Escherichia coli
<b>Protein length</b>	Full length protein
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Amino acids</b>	21 to 677

### Specifications

Our **Abpromise guarantee** covers the use of **ab93930** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	SDS-PAGE
<b>Form</b>	Liquid

### Preparation and Storage


<b>Stability and Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
------------------------------	---

### General Info

<b>Function</b>	Secretogranin-1 is a neuroendocrine secretory granule protein, which may be the precursor for other biologically active peptides.
<b>Tissue specificity</b>	Expressed in the adrenal medulla, and in pheochromocytoma. Not expressed in liver.

<b>Sequence similarities</b>	Belongs to the chromogranin/secretogranin protein family.
<b>Post-translational modifications</b>	Extensively processed by limited proteolysis at conserved basic residues. Alternative processing are seen in different tissues.
<b>Cellular localization</b>	Secreted. Neuroendocrine and endocrine secretory granules.

## Images



15% SDS-PAGE showing ab93930 (3 µg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors