

Product datasheet

Recombinant Human CKS1 protein (His tag) ab177620

[1 Image](#)

Description

Product name	Recombinant Human CKS1 protein (His tag)
Purity	> 95 % SDS-PAGE. purified by using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>P61024</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMSHKQIY YSDKYDDEEF EYRHVMLPKD IAKLVPKTHL MSESEWRNLG VQQSQGWVHY MIHEPEPHIL LFRRPLPKKP KK
Predicted molecular weight	12 kDa including tags
Amino acids	1 to 79
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab177620** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------------------------	--

pH: 8.00


Constituents: 0.32% Tris HCl, 0.88% Sodium chloride, 20% Glycerol (glycerin, glycerine), 0.02% DTT

General Info

Function Binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function.

Sequence similarities Belongs to the CKS family.

Images



15% SDS-PAGE analysis of ab177620 (3µg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors