

Recombinant Human CPEB3 protein ab153194

1 Image

Description	
Product name	Recombinant Human CPEB3 protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<div>MQDDLMDKSKTQPQPQQQQRQQQPQPESSVSEAPS</div> <div>TPLSSETPKPEEN</div> <div>SAVPALSPAAAPPAPNGPDKMQMESPLLPGLSFHQPPQ</div> <div>QPPPPQEPAAPG</div> <div>ASLSPSFGSTWSTGTTNAVEDSFFQGITPVNGTMLFQNFP</div> <div>HHVNPVFGGT</div> <div>FSPQIGLAQTQHHQQPPPPAPAPQPAQPAQPPQAQPPQ</div> <div>QRRSPASPSQAP</div> <div>YAQRSAAAAYGHQPIMTSKPSSSSAVAAAAAAAAAASSAS</div> <div>SSWNTHQSVNA</div> <div>AWSAPSNPWGGLQAGRDPRRAVGVGVGVGVPSPLN</div> <div>PISPLKKPFSSNV</div> <div>IAPPKFRAAPLTSKSWMEDNAFWTDNGNNLLPFQDRSR</div> <div>PYDTFNLHSLE</div> <div>NSLMDMIRTDHEPLKGKHYPSPGPPMSFADIMWRNHFAG</div> <div>RMGINFHHPGT</div> <div>DNIMALNNAFLDDSHGDQALSSGLSSPTRCQNGERVERY</div> <div>SRKV FVGGLPP</div> <div>DIDEDEITASFRRFGPLVVDWPHKAESKSYFPPKGYAFL</div> <div>FQEESSVQAL</div> <div>IDACLEEDGKLYLCVSSPIKDKPVQIRPWNLSDSDFVMDG</div> <div>SQPLDPRKTI</div> <div>FVGGVPRPLRAVELAMIMDRLYGGVCYAGIDTDPCLKYPK</div> <div>GAGRVAFSNQ</div> <div>QSYAAISARFVQLQHNDIDKRVEVKPYVLDDQMCDECQG</div> <div>TRCGGKFAPF</div> <div>FCANVTCLQYYCEYCWASIHSRAGREFHKPLVKEGGDRP</div> <div>RHVPFRWS</div>

<b>Amino acids</b>	1 to 697
<b>Tags</b>	GST tag N-Terminus

## Specifications

---

Our **Abpromise guarantee** covers the use of **ab153194** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** Western blot

ELISA

**Form** Liquid

**Additional notes**

## Preparation and Storage

---

**Stability and Storage** Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl


## General Info

---

**Sequence similarities** Belongs to the RRM CPEB family.  
Contains 2 RRM (RNA recognition motif) domains.

## Images

---



ab153194 on a 12.5% SDS-PAGE stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors