

Product datasheet

Recombinant Human DcR3 protein ab157279

1 Image

Description

Product name	Recombinant Human DcR3 protein	
Biological activity	ab157279 binds Human LIGHT. Does not bind Human FasL.	
Purity	> 95 % SDS-PAGE.	
Endotoxin level	< 0.100 Eu/μg	
Expression system	HEK 293 cells	
Accession	<u>O95407</u>	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MRALEGPGLS LLCLVLALPA LLPVPAVRGV AETPTYPW RD AETGERLVCA QCPPGTFVQR PCRRDSPTTC GPCPPRHYTQ FWNLERCRY CNVLCGEREE EARACHATHN RACRCRTGFF AHAGFCLEHA SCPPGAGVIA PGTPSQNTQC QPCPPGTFSA SSSSSEQCQP HRNCTALGLA LNVPGSSSHD TLCTSCTGFP LSTRVPGAEE CERAVIDFVA FQDISIKRLQ RLLQALEAPE GWGPTP	
Predicted molecular weight	58 kDa including tags	
Amino acids	1 to 246	

Specifications

Our **Abpromise guarantee** covers the use of **ab157279** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Functional Studies SDS-PAGE
Form	Lyophilized
Additional notes	ab157279 binds Human LIGHT. Does not bind Human FasL.

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at -20°C.

Constituent: 99% PBS

Reconstitution

Reconstitute with 50µl sterile water to a concentration of 1 mg/ml. After reconstitution, prepare aliquots and store at -20°C. Avoid freeze/thaw cycles.

General Info

Function

Decoy receptor for the cytotoxic ligands TNFS14/LIGHT and TNFSF6/FASL. Protects against apoptosis.

Tissue specificity

Detected in fetal lung, brain and liver. Detected in adult stomach, spinal cord, lymph node, trachea, spleen, colon and lung. Highly expressed in several primary tumors from colon, stomach, rectum, esophagus and in SW480 colon carcinoma cells.


Sequence similarities

Contains 4 TNFR-Cys repeats.

Cellular localization

Secreted.

Images



ELISA analysis of ab157279.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors